



Project TA01011423

Hamry 3567, 580 01 Havl. Brod, Czech Republic e-mail: 1-cube @ 1-cube.com web site: www.1-cube.com

Laboratory Hazemeter/Colorimeter COLORTURB- MZN2012

Handy design verified by a daily laboratory practice, extremely easy operation, modern optoelectronic components, microprocessor control – all of this makes the hazemeter useful for everyday use in any place, where the clarity/haze and color of liquid is to be handled regularly. Extremely long lifetime (up to 10 years, 200 000 measurements)

Routine laboratory check of products and intermediates



food industry water supply - drinking water treatment - wastewater processing pharmaceutical and cosmetic industry paints and inks production petrochemical industry photography and films processing textile industry energetics electrotechnical industry

Immerse liquid technique, automatic rotating of a sample and statistical processing of measured signals easily allows measurement of samples not only in standard cuvettes, but directly in closed bottles of various diameter as well.

- Direct quality check of products sealed in bottles
- Shelf live prediction and control



- beverage industry brewery
 viniculture
 inice and fresh drinks proceed.
 - juice and fresh drinks production

The gathering information power may be amplified many times after connecting to PC. The special controlling SW together with various accessories and adapters move function of this hazemeter far beyond the common ones and allows performing wide range of tasks.

Wireless adapter (Bluetooth serial adapter from Socket) enables the customer to control several hazemeters from one PC.


- Measurement of small volume samples in chemical tubes
- Kinetic measurement of haze development
- Turbidimetric titrations
- Wireless bluetooth connection

Principle of operation:





The hazemeter/color MZN based on the light scattering technique is used to detect the haze particles present in monitored liquid. The light coming from the source is focused by the system of apertures and lenses into a narrow beam entering a sample. A part of the beam is scattered around by haze



particles being in a liquid. There are two optical detectors set in the hazemeter MZN, that measure simultaneously the light scattered in two angles -90° and 15° (so called N - nephelometric and F- forward). The third optical detector evaluates the light intensity changes of the beam transmitted through a sample (0° T-transmission). The both scattering signals N, F are rationed to the transmission one (in the ratio mode). This provides long time calibration stability and compensates the light attenuation caused by possible coloration of a sample (colour compensation).

The fourth detector R – Reference monitors optical source stability. Ratio Transmision to reference is proportional to light attenuation in the sample (Color).

Calibration relationships for haze based on formazin haze standard are stored in the memory of the hazemeter. Measured values for all the channels (H90 and H12 and Absorbance) are then displayed in chosen haze units on the LCD display.

Dual-angle hazemeter offers qualitatively entirely new possibilities compared to common single-angle nephelometers or turbidimeters. Intensity of light scattered at specific angle depends not only on the concentration of haze particles but also on their size, therefore it is possible to gather information regarding the particle size distribution.

Dual-angle haze and light attenuation (Color) are measured independently on three optical wavelengths (Red, Green, Blue).

	Characteristics	Technical	parameters
		optical set-up	dual-angle hazemeter with transmission
٠	easy operation, just one touch button operation		compensation
٠	immediate readiness to measure	measuring angles	90° , 15° , transmission (0°)
•	measurement in cuvettes, chemical tubes and	measuring range for H90	0 - 250 EBC (1000 EBC on demand)
	bottles	H12	0 - 80 EBC
٠	automatic bottle rotating	Abs	0 – 5 Abs
٠	short measuring time	units	EBC, NTU, ppm, (other on request)
٠	digital signal processing	ranging	fully automatic
٠	long-term stability	resolution	1% on LCD display
٠	easy recalibration		0.1 % in digital output
٠	fully automatic range switching	accuracy	2% of measured value 0,01 up to 100EBC
٠	simple and practical design		5% of measured value up to 250 EBC
٠	stainless steel upper casing	reproducibility	1 % of measured value
٠	clearly visible LCD display	immersion liquid	distilled water
•	PC connectivity	measuring time	7 – 12 s
•	sample thermostating possibility		
		light source	LED RGB
sep	arate accessories:	- colour spectrum	Red (650nm), Green (530nm),
-	adapter for chemical tubes		Blue (450nm)
-	adapter for chemical tubes with magnetic stirrer	- lifetime	100 000 measurements
-	flow-through measurement cells thermo-electrical thermostat, thermometer	sample volume	from 150 ml in cuvette
-	dosing pumps		from 5 ml in chemical tube
-	PC software	di sital suturet	bottles 0,3 – 2,0litre (diameter 100mm) RS 232, Bluetooth
-	- kinetics of haze development	digital output	,
	- turbidimetric titrations	display	16 characters two lines illuminated LCD display
	user-defined calibrationsother special features on request	voltage, power input	230 V/50 Hz, 50 VA
		dimensions, weight	450x350x190 mm, 9 kg
		chamber diameter	100 mm