



# 1-CUBE

Hamry 3567, 58001 Havlíčkův Brod, Czech Republic.

tel + 420 569 433 620

1-cube@1-cube.com

www.1-cube.com

## DECARBONIZER – TYPE OP

### Application:

- **Efficient sample preparation:** Rapid removal of CO<sub>2</sub> from the sample for subsequent chemical analysis.
- **Real-time fermentation monitoring:** Immediate insight into the fermentation process in CCTs, maturation tanks, and other tanks.

### Technical description:

Most chemical analyses of beer and carbonated beverages require prior removal of CO<sub>2</sub>. The OP-type decarbonizer prepares a 150 ml sample for analysis within one minute. In combination with an aseptic sampling valve (type AOV) and a saccharometer, it forms a complete system for fast and accurate determination of the degree of fermentation directly from process tanks.

The device is designed for very rapid decarbonization of beer and carbonated beverage samples based on the principle of mechanical stirring. The assembly consists of a transparent plexiglass container and a removable lid with an integrated stirrer. The control electronics allow precise adjustment of the stirring time, making it possible to define the required level of CO<sub>2</sub> removal.



**Technical data:**

<b>Parameter</b>	<b>Range / Value</b>
Maximum volume of decarbonized sample	150 ml
Maximum CO <sub>2</sub> content after 60 s of decarbonization	0.05%
Power supply	230 V / 50 Hz, 2.5 A
Dimensions (L × W × H)	320 × 320 × 320 mm
Weight	4 kg

**Scope of Delivery:**

- Decarbonizer type OP, User Manual

**Accessories (optional):**

- glass graduated test tube
- brewing hydrometers

**Advantages & Benefits:**

- Very fast removal of CO<sub>2</sub> from the beverage sample prior to analysis.
- Compared to conventional shakers, CO<sub>2</sub> removal from the sample is an order of magnitude faster.

**FAQ:****How to quickly determine the degree of fermentation using the OP decarbonizer?**

The measurement procedure is fast and requires only three simple steps:

1. **Sample collection:** Take approximately 150 ml of beer from the cylindroconical tank (CCT).
2. **Decarbonization:** Place the sample into the OP device. Within 1 minute, complete removal of CO<sub>2</sub> from the sample is achieved, which is essential for measurement accuracy.
3. **Measurement:** Pour the degassed beer into the graduated measuring tube (integrated directly in the decarbonizer) and immediately read the current degree of fermentation using a brewing hydrometer (saccharometer).